



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Cadmium(Cd) Monoclonal Antibody

Product name: Anti-Cadmium Monoclonal Antibody

Code: LD-ZJS-1001-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 4-5ppb

CollaurumSensitivity: 25ppb

Protein Concentration: 8mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Lead(Pb) monoclonal antibody, anti-Mercury(Hg) monoclonal antibody, anti-Methyl Mercury Chloride(CH₃HgCl) monoclonal antibody, anti-Copper(Cu) monoclonal antibody, anti-Chromium(Cr) monoclonal antibody.