



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Chromium(Cr) Monoclonal Antibody

Product name: Anti-Chromium(Cr) Monoclonal Antibody

Code: LD-ZJS-1006-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 18-20ppb

CollaurumSensitivity: 1000ppb

Protein Concentration: 8mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Cadmium monoclonal antibody, anti-Lead(Pb) monoclonal antibody, anti-Mercury(Hg) monoclonal antibody, anti-Methyl Mercury Chloride monoclonal antibody, anti-Copper(Cu) monoclonal antibody.