



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Anti-Mercury(Hg) Monoclonal Antibody**

**Product name:** Anti-Mercury(Hg) Monoclonal Antibody

**Code:** LD-ZJS-1003-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 4-5ppb

**CollaurumSensitivity:** 250ppb

**Protein Concentration:** 8mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Cadmium monoclonal antibody, anti-Lead(Pb) monoclonal antibody, anti-Methyl Mercury Chloride(CH<sub>3</sub>HgCl) monoclonal antibody, anti-Copper(Cu) monoclonal antibody, anti-Chromium(Cr) monoclonal antibody.