



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Furaltadone Metabolite Monoclonal Antibody

Product name: Anti-Furaltadone Metabolite Monoclonal Antibody

Code: LD-XY-1003-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 1ppb

ELISA Sensitivity: 0.1ppb

CollaurumSensitivity: 10ppb

Protein Concentration: 8mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Tolfenamicacid monoclonal antibody, anti-Carprofen monoclonal antibody, anti-Hydrocortisone monoclonal antibody, anti-Metamizole Sodium Tablets monoclonal antibody, anti-Clonazepam monoclonal antibody, anti-Diazepam monoclonal antibody, anti-Antipyrine monoclonal antibody, anti-Indometacin Tablets monoclonal antibody.

.

.