



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Anti-Metamizole Sodium Tablets Monoclonal Antibody**

**Product name:** Anti-Metamizole Sodium Tablets Monoclonal Antibody

**Code:** LD-XY-1005-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 2ppb

**ELISA Sensitivity:** 0.2ppb

**CollaurumSensitivity:** 18ppb

**Protein Concentration:** 7mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Tolfenamicacid monoclonal antibody, anti-Carprofen monoclonal antibody, anti-Furaltadone monoclonal antibody, anti-Hydrocortisone monoclonal antibody, anti-Clonazepam monoclonal antibody, anti-Diazepam monoclonal antibody, anti-Antipyrine monoclonal antibody, anti-Indometacin Tablets monoclonal antibody.

.

.