



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Anti-Clorprenaline Monoclonal Antibody(CLP)**

**Product name:** Anti-Clorprenaline Monoclonal Antibody

**Code:** LD-SY-1007-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 5ppb

**ELISA Sensitivity:** 0.5ppb

**CollaurumSensitivity:** 20ppb

**Protein Concentration:** 5mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Ractopamine monoclonal antibody, anti-Salbutamol monoclonal antibody, anti-Chlorpromazine monoclonal antibody, anti-Cimaterol monoclonal antibody, anti-Cyproheptadine monoclonal antibody.