



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Phenylethylamine Monoclonal Antibody(PEA)

Product name: Anti-Phenylethylamine Monoclonal Antibody

Code: LD-SY-1003-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.5ppb

ELISA Sensitivity: 0.1ppb

CollaurumSensitivity: 20ppb

Protein Concentration: 5mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Ractopamine monoclonal antibody, anti-Salbutamol monoclonal antibody, anti-Chlorpromazine monoclonal antibody, anti-Cimaterol monoclonal antibody, anti-Cyproheptadine monoclonal antibody.