



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-25-OH-Vitamin D3 Monoclonal Antibody

Product name: Anti-25-OH-Vitamin D3 Monoclonal Antibody

Code: LD-SLHX-1005-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.4ppb

CollaurumSensitivity: 250ppb

Protein Concentration: 6mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Folate monoclonal antibody, anti-Vitamin B12 monoclonal antibody.