



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

## **Anti-Aflatoxin B1 Monoclonal Antibody**

**Product name:** Anti-Aflatoxin B1 Monoclonal Antibody

**Code:** LD-DS-1001-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 1ppb

**CollaurationSensitivity:** 10ppb

**Protein Concentration:** 5mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-T-2 toxin monoclonal antibody, anti-Aflatoxin M1 monoclonal antibody, anti-Total aflatoxin monoclonal antibody, anti-Ochratoxin A monoclonal antibody, anti-Zearalenone monoclonal antibody, anti-Fumonisin monoclonal antibody, anti-Deoxynivalenol monoclonal antibody, anti-Citrinin monoclonal antibody, anti-Sterigmatocystin monoclonal antibody, anti-Alternariol monoclonal antibody, anti-Ricin monoclonal antibody, anti-Microcystin LR monoclonal antibody, anti-Diphtheria toxin monoclonal antibody, anti-Cholera toxin monoclonal antibody, anti-Clostridium botulinum toxinA/B/D/E monoclonal antibody, anti-E. coli heat-labile enterotoxin A、B monoclonal antibody.