



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Aflatoxin M1 Monoclonal Antibody

Product name: Anti-Aflatoxin M1 Monoclonal Antibody

Code: LD-DS-1002-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 1ppb

CollaurumSensitivity: 0.5ppb

Protein Concentration: 8mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-T-2 toxin monoclonal antibody, anti-Aflatoxin B1 monoclonal antibody, anti-Total aflatoxin monoclonal antibody, anti-Ochratoxin A monoclonal antibody, anti-Zearalenone monoclonal antibody, anti-Fumonisin monoclonal antibody, anti-Deoxynivalenol monoclonal antibody, anti-Citrinin monoclonal antibody, anti-Sterigmatocystin monoclonal antibody, anti-Alternariol monoclonal antibody, anti-Ricin monoclonal antibody, anti-Microcystin LR monoclonal antibody, anti-Diphtheria toxin monoclonal antibody, anti-Cholera toxin monoclonal antibody, anti-Clostridium botulinum toxinA/B/D/E monoclonal antibody, anti-E. coli heat-labile enterotoxin A、B monoclonal antibody.