



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Human Chorionic Gonadotropin Monoclonal Antibody

Product name: Anti-Human chorionic gonadotropin Monoclonal Antibody

Code: LD-JS-1012-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.2ppb

CollaurationSensitivity: 10ppb

Protein Concentration: 6mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Erythropoietin monoclonal antibody, anti-Thyroxin monoclonal antibody, anti-Calcitonin monoclonal antibody, anti-Dexamethasone monoclonal antibody, anti-Cortisol monoclonal antibody, anti-Estrone monoclonal antibody, anti-Progesterone monoclonal antibody, anti-Estrone monoclonal antibody, anti-Estradiol monoclonal antibody, anti-Growth hormone monoclonal antibody, anti-Prostaglandin monoclonal antibody, anti-uteinizing hormone/Prolactin monoclonal antibody, anti-Bisphenol A monoclonal antibody, anti-Testosterone monoclonal antibody.