



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

## **Anti-Prostaglandin Monoclonal Antibody**

**Product name:** Anti-Prostaglandin Monoclonal Antibody

**Code:** LD-JS-1005-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 1-2ppb

**CollaurumSensitivity:** 100ppb

**Protein Concentration:** 6mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Erythropoetin monoclonal antibody, anti-Thyroxin monoclonal antibody, anti-Calcitonin monoclonal antibody, anti-Dexamethasone monoclonal antibody, anti-Cortisol monoclonal antibody, anti-Estrone monoclonal antibody, anti-Progesterone monoclonal antibody, anti-Estrone monoclonal antibody, anti-Estradiol monoclonal antibody, anti-Growth hormone monoclonal antibody, anti-Human chorionic gonadotropin monoclonal antibody, anti-uteinizing hormone/Prolactin monoclonal antibody, anti-Bisphenol A monoclonal antibody, anti-Testosterone monoclonal antibody.