



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Uteinizing Hormone、 Prolactin Monoclonal Antibody

Product name: Anti-Uteinizing Hormone、 Prolactin Monoclonal Antibody

Code: LD-JS-1008-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.5ppb

CollaurumSensitivity: 25ppb

Protein Concentration: 6mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Dexamethasone monoclonal antibody, anti-Calcitonin monoclonal antibody, anti-Cortisol monoclonal antibody, anti-Testosterone monoclonal antibody, anti-Prostaglandin monoclonal antibody, anti-Erythropoetin monoclonal antibody, anti-Progesterone monoclonal antibody, anti-Estrone monoclonal antibody, anti-Estradiol monoclonal antibody, anti-Growth hormone monoclonal antibody, anti-Human chorionic gonadotropin monoclonal antibody, anti-Diethylstilbestrol/Hexestrol monoclonal antibody, anti-Bisphenol A monoclonal antibody, anti-Thyroxin(e) monoclonal antibody.