



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-2-Nitrofurantoin metabolite Monoclonal Antibody

Product name: Anti-2-Nitrofurantoin metabolite Monoclonal Antibody

Code: LD-KSS-1060-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.6ppb

ELISA Sensitivity: 0.1ppb

CollaurumSensitivity: 5ppb

Protein Concentration: 6mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-2-Nitrofurazone metabolite monoclonal antibody, anti-2-Furaltadone metabolite monoclonal antibody, anti-2-Furazolidone metabolite monoclonal antibody.

.

.