



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti-Ampicillin Monoclonal Antibody

Product name: Anti-Ampicillin Monoclonal Antibody

Code: LD-KSS-1046-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 20ppb

ELISA Sensitivity: 2ppb

CollaurumSensitivity: 1000ppb

Protein Concentration: 8mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Ceftiofur monoclonal antibody, anti-Cefquinome monoclonal antibody, anti-Cephalosporin monoclonal antibody, anti-Cephalexin monoclonal antibody, anti-Hygromycin B monoclonal antibody, anti-Bacitracin monoclonal antibody, anti-Polymyxin B monoclonal antibody, anti-Ampicillin monoclonal antibody, anti-Cyclosporin A monoclonal antibody, anti-metoclopramide monoclonal antibody, anti-Trimethoprim monoclonal antibody, anti-Polymyxin E monoclonal antibody.

.

.