



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Anti- Ampicillin Trihydrate Monoclonal Antibody**

**Product name:** Anti- Ampicillin Trihydrate Monoclonal Antibody

**Code:** LD-KSS-1054-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 2ppb

**ELISA Sensitivity:** 0.2ppb

**CollaurumSensitivity:** 10ppb

**Protein Concentration:** 3.8mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Ampicillin monoclonal antibody, anti-Ceftiofur monoclonal antibody, anti-Cefquinome monoclonal antibody, anti-Cephalosporin monoclonal antibody, anti-Cephalexin monoclonal antibody, anti-Hygomycin B monoclonal antibody, anti- Bacitracin monoclonal antibody, anti- Polymyxin B monoclonal antibody, anti-Cyclosporin A monoclonal antibody, anti-metoclopramide monoclonal antibody, anti-Trimethoprim monoclonal antibody, anti-Polymyxin E monoclonal antibody.

.

.