



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti- Cyclosporin A Monoclonal Antibody

Product name: Anti- Cyclosporin A Monoclonal Antibody

Code: LD-KSS-1055-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 2ppb

ELISA Sensitivity: 0.2ppb

CollaurumSensitivity: 10ppb

Protein Concentration: 7.3mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Ampicillin monoclonal antibody, anti-Ceftiofur monoclonal antibody, anti-Cefquinome monoclonal antibody, anti-Cephalosporin monoclonal antibody, anti-Cephalexin monoclonal antibody, anti-Hygomycin B monoclonal antibody, anti- Bacitracin monoclonal antibody, anti- Polymyxin B monoclonal antibody, anti- Ampicillin Trihydrate monoclonal antibody, anti-metoclopramide monoclonal antibody, anti-Trimethoprim monoclonal antibody, anti-Polymyxin E monoclonal antibody.

.

.