



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

Anti- Amikacin Monoclonal Antibody

Product name: Anti- Amikacin Monoclonal Antibody

Code: LD-KSS-1012-2

Source: BALB/C mouse

Product Type: IgG2a

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:400000

IC50: 0.5ppb

ELISA Sensitivity: 0.1ppb

CollaurumSensitivity: 10ppb

Protein Concentration: 6mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 5 years, kept below 70 °C

Related products: anti-Kanamycin monoclonal antibody, anti-Gentamicin monoclonal antibody, anti- Streptomycin monoclonal antibody, anti- Apramycin monoclonal antibody, anti-Spectinomycin monoclonal antibody, anti- Tobramycin monoclonal antibody, anti-Neomycin monoclonal antibody.

.

.