



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Anti- Apramycin Monoclonal Antibody**

**Product name:** Anti- Apramycin Monoclonal Antibody

**Code:** LD-KSS-1014-2

**Source:** BALB/C mouse

**Product Type:** IgG2a

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:400000

**IC50:** 5ppb

**ELISA Sensitivity:** 0.5ppb

**CollaurumSensitivity:** 100ppb

**Protein Concentration:** 5mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 5 years, kept below 70 °C

**Related products:** anti-Kanamycin monoclonal antibody, anti-Gentamicin monoclonal antibody, anti- Amikacin monoclonal antibody, anti- Streptomycin monoclonal antibody, anti- Spectinomycin monoclonal antibody, anti- Tobramycin monoclonal antibody, anti-Neomycin monoclonal antibody.

.

.