



**Shandong Lvdu Bio-Sciences&Technology Co., Ltd.**

---

**Clostridium Botulinum ToxinA/B/D/E Antigen**

**Product name:** Clostridium Botulinum ToxinA/B/D/E Antigen

**Code:** LD-DS-1016-1

**Protein carrier:** BALB/C mouse

**Product type:** BSA/OVA

**Physical properties:** Liquid or solid, avoid repeated freezing and thawing

**ELISA titer:** 0.005M PBS, > 1:30000

**Collaurum titer:** > 1:40, coupling ratio > 1:60

**Protein concentration:** 8.0mg/ml

**Buffer:** 0.01M PBS PH7.4

**Storage:** Below- 20 °C

**Validity:** 2 years, avoid repeated freezing and thawing