



Shandong Lvdu Bio-Sciences&Technology Co., Ltd.

2,4,5-Trichlorophenoxyacetic Antigen

Product name: 2,4,5-Trichlorophenoxyacetic Antigen

Code: LD-NY-1016-1

Protein carrier: BALB/C mouse

Product type: BSA/OVA

Physical properties: Liquid or solid, avoid repeated freezing and thawing

ELISA titer: 0.005M PBS, > 1:30000

Collaurum titer: > 1:40, coupling ratio > 1:60

Protein concentration: 8.0mg/ml

Buffer: 0.01M PBS PH7.4

Storage: Below- 20 °C

Validity: 2 years, avoid repeated freezing and thawing